

Clear benefits

Diagnose faster and more accurately

- Match look & feel to your diagnostic needs
- One-touch filter switching
- Get the image you need more quickly

Reduce the risk of misdiagnosis

- See clinical details better
- Improve certainty by comparing views
- No artifacts to cloud the picture

Work the way you want to

- Choose from up to six look & feel options
- Create your own filter settings
- Intuitive software and real-time filtering

Always count on the same experience

- Rendering algorithms deliver same great results across devices
- Works with CS 8100, CS 8100 3D, CS 9000 Family and CS 9300 Family
- Same quality for panoramic and cephalometric images

Let's redefine expertise



CS 9000/CS 9300

CS 8100

CS 8100 3D



CS ADAPT IS AVAILABLE FOR ALL CARESTREAM
DENTAL EXTRAORAL IMAGING SYSTEMS*

* Not available for CS 8000C


HOW TO UPGRADE

Call XXX XXX XX or contact your authorized dealer
to order your software upgrade today.

Discover more at carestreamdental.com



A whole new level of choice
in extraoral image processing



THE WORLD'S ONLY
CUSTOMIZABLE LOOK & FEEL

Upgrade to more individualized diagnostics

Imaging that works the way you do

You don't work in exactly the same way as other dentists, so why should your extraoral imaging system give you the exact same images it gives everyone else? With the CS Adapt image processing module, it doesn't any more.

This major upgrade to our Dental Imaging Software lets you decide how you want your images to appear, by choosing from up to six different look & feel options for your panoramic and cephalometric images.

In addition to letting you choose the look & feel of your images, the CS Adapt upgrade features state-of-the-art rendering algorithms that dramatically improve image quality across the Carestream Dental range.

So whether you're working with 2D panoramic or cephalometric images, CS Adapt will let you work with greater confidence and less risk of misdiagnosis.



CS ADAPT filter library



Panoramic:

ADAPT ORIGINAL



The look & feel you know from your current system

ADAPT CONTRAST



Optimizes contrast for dense structures like cortical bone

ADAPT T-MAT



Recreates a familiar X-ray film quality

ADAPT DYNAMIC



Our highest contrast filter

ADAPT SHARP



High contrast distinguishes between dentine and alveolar bone

ADAPT SMOOTH



Brings out soft tissue with medium contrast

Software that makes a clear difference



Upgrading your Carestream Dental Imaging Software with CS Adapt module gives you:

- ✓ Six panoramic look & feel filter presets
- ✓ Four cephalometric look and feel filter presets
- ✓ Ability to define your own look & feel

- ✓ Intuitive library-style browsing and selection
- ✓ State-of-the-art image quality
- ✓ No more processing artifacts or "dark halos"
- ✓ Same flexibility and quality across devices

CS ADAPT filter library

Cephalometric:

ADAPT ORIGINAL



The look & feel you know from your current system

ADAPT SHARP



Enhances visibility with high sharpness and contrast

ADAPT T-MAT



Recreates a familiar X-ray film quality

ADAPT SMOOTH



Enhances visibility with sharpness and medium contrast